



CATALOGUE – August 2016 update

**Walkers
Walking Frames
Walking Aids
Wheelchairs
Scooters
Bedroom
Powerlift Chairs
Bathroom
Aids to Daily Living**



*** Also Hire, Service & Repair of our Products**

**ADDRESS: 473 Ballarat Road
Sunshine Vic 3020**

PHONE: (03) 9312 2888

FAX: (03) 9312 2999

EMAIL: enquiries@wheelchairman.com.au

WEB: www.wheelchairman.com.au

**See map on
back cover**

WALKERS - STANDARD SIZES

A-FRAME WALKER. ALUMINIUM. CABLE BRAKES. FIXED, REMOVABLE, CURVED, BACKREST WM1003WCM

Comfortable padded seat
Semi-curved backrest
Removable storage basket and bag
Adjustable handle height
Folds into compact unit for easy storage
Seat Height: 56 cm
Seat Width: 45 cm
Seat Depth: 37 cm
Total Width: 64 cm
Wheel Size: 17.5 cm (7")
Handle Height: 84 - 98 cm
Max. Weight Capacity: 130 kg
Weight: 7.9 kg
Colour: Dark green



A-FRAME WALKER. ALUMINIUM. CABLE BRAKES. SWIVEL BACKREST WM1003ND

Seat Height: 56 cm
Seat Width: 46 cm
Seat Depth: 33 cm
Total Width: 62 cm
Wheel Size: 19 cm (7 1/2")
Handle Height: 79 - 94 cm
Max. Weight Capacity: 130 kg
Weight: 7.6 kg
Colour: Blue / Burgundy / Titanium



A-FRAME WALKER. ALUMINIUM. CABLE BRAKES. FIXED / SWIVEL BACKREST WM1003P

Seat Height: 56 cm
Seat Width: 45 cm
Seat Depth: 31 cm
Total Width: 62 cm
Wheel Size: 15 cm (6")
Handle Height: 82 - 93 cm
Max. Weight Capacity: 100 kg
Weight: 6.3 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

**A-FRAME WALKER. ALUMINIUM.
WEIGHT-ACTIVATED BRAKES.
WM1002**

Comfortable padded seat
Removable storage basket
Powder coated aluminium
Folds into compact unit for easy storage
Seat Height: 56 cm
Seat Width: 45 cm
Seat Depth: 30 cm
Total Width: 63 cm
Wheel Size: 17.5 cm (7")
Handle Height: 84-94 cm
Max. Weight Capacity: 120 kg
Weight: 7.4 kg



**A-FRAME WALKER. STEEL (RIVA). CABLE BRAKES
WM1004**

4 Wheel steel walker
Adjustable cable brakes
Basket and tray included
Seat Height: 60 cm
Seat Width: 43 cm
Seat Depth: 22 cm
Total Width: 60 cm
Wheel Size: 200 x 50 mm
Handle Height: 78 – 95 cm
Max. Weight Capacity: 100 kg
Weight: 12 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

SMALLER SIZES

**A-FRAME WALKER. ALUMINIUM. CABLE BRAKES.
FIXED, REMOVABLE, CURVED, BACKREST
WM1012WCM**

Seat Height: 49 cm
Seat Width: 39.5 cm
Seat Depth: 37 cm
Total Width: 59 cm
Wheel Size: 15 cm (6")
Handle Height: 76 – 89 cm
Max. Weight Capacity: 100 kg
Weight: 7.1 kg



**A-FRAME WALKER. ALUMINIUM. CABLE BRAKES.
FIXED, REMOVABLE, CURVED, BACKREST
WM1012ND**

Seat Height: 47 cm
Seat Width: 46 cm
Seat Depth: 32 cm
Total Width: 61 cm
Wheel Size: 15 cm (6")
Handle Height: 71 – 82 cm
Max. Weight Capacity: 130 kg
Weight: 6.25 kg



**A-FRAME WALKER. ALUMINIUM. CABLE BRAKES.
SMALL FRAME. MEDIUM-WIDE.
FIXED / SWIVEL BACKREST
WM1012MW**

Seat Height: 50 cm
Seat Width: 45 cm
Seat Depth: 32 cm
Total Width: 62 cm
Wheel Size: 19 cm (7 1/2")
Handle Height: 75 – 85 cm
Max. Weight Capacity: 130 kg
Weight: 7.4 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

**A-FRAME WALKER. ALUMINIUM. CABLE BRAKES.
FIXED / SWIVEL BACKREST
WM1012SW**

Seat Height: 47 cm
Seat Width: 40 cm
Seat Depth: 30cm
Total Width: 57 cm
Wheel Size: 15 cm (6")
Handle Height: 74 – 84 cm
Max. Weight Capacity: 100 kg
Weight: 6.2 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

WALKERS - BARIATRIC MODELS

**BARIATRIC A-FRAME WALKER. ALUMINIUM.
CABLE BRAKES. FIXED, REMOVABLE, CURVED,
BACKREST. BAG AND BASKET UNDER SEAT
WM1003LWCM**

Seat Height: 58 cm
Seat Width: 52 cm
Seat Depth: 36 cm
Total Width: 72 cm
Wheel Size: 17.5 cm (7")
Handle Height: 85 – 98.5 cm
Max. Weight Capacity: 150 kg
Weight: 8.1 kg



**BARIATRIC A-FRAME WALKER. ALUMINIUM.
CABLE BRAKES. FIXED STRAIGHT BACKREST.
BASKET IN FRONT
WM1003L**

Seat Height: 58 cm
Seat Width: 52 cm
Seat Depth: 32 cm
Total Width: 63 cm
Wheel Size: 20 cm (8")
Handle Height: 88 – 99 cm
Max. Weight Capacity: 150 kg
Weight: 10 kg



**BARIATRIC A-FRAME WALKER. ALUMINIUM. CABLE BRAKES. FIXED STRAIGHT
BACKREST. BASKET IN FRONT
WM1003XL – WIDER MODEL**

Seat Height: 56 cm
Seat Width: 59 cm
Seat Depth: 35 cm
Total Width: 71 cm
Wheel Size: 19 cm (7 1/2")
Handle Height: 77 – 93 cm
Max. Weight Capacity: 200 kg
Weight: 11.6 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

**BARIATRIC A-FRAME WALKER. ALUMINIUM. CABLE
BRAKES. FIXED, REMOVABLE, CURVED, BACKREST.
BAG & BASKET UNDER SEAT
WM1012WWCM – SMALL**

Seat Height: 49 cm
Seat Width: 52 cm
Seat Depth: 36 cm
Total Width: 72 cm
Wheel Size: 17.5 cm (7")
Handle Height: 75 – 87 cm
Max. Weight Capacity: 150 kg
Weight: 8 kg



**BARIATRIC A-FRAME WALKER. ALUMINIUM. CABLE BRAKES.
FIXED STRAIGHT BACKREST. BASKET IN FRONT
WM1012W – SMALL**

Seat Height: 50 cm
Seat Width: 52 cm
Seat Depth: 32 cm
Total Width: 63 cm
Wheel Size: 20 cm (8")
Handle Height: 75 – 86 cm
Max. Weight Capacity: 150 kg
Weight: 9.9 kg



**OXYGEN CYLINDER HOLDER TO FIT 4 WHEEL WALKERS
VARIOUS SIZES
WM1021**



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

**TRI-WALKER. ALUMINIUM. CABLE BRAKES INCLUDING
BAG & BASKET
WM1010**

3 Wheel walker ideal for small areas or indoor use
Total Width: 64 cm
Wheel Size: 20 x 5 cm solid tyre
Handle Height: 81 – 94 cm
Max. Weight Capacity: 100 kg
Weight: 7 kg



EASY WALKERS

**EASY WALKER. FOLDING WITH FOREARM SUPPORTS.
FRONT WHEELS, SKIS / FERRULES
WM1018F**

2 x 3" Wheel walking frame with rear skis / ferrules
Total Width: 60 cm
Total Length: 66 cm
Handle Height: 95 – 125 cm
Max. Weight Capacity: 100 kg.
Weight: 7.1 kg



WM1018NF Also available in NON-FOLDING

**EASY WALKER. FOLDING WITH FOREARM SUPPORTS.
CABLE BRAKES. 4 WHEELS.
WM1018BS**

Total Width: 61 cm
Total Length: 69 cm
Wheel Size: 200 x 50 mm
Handle Height: 97 – 114 cm
Max. Weight Capacity: 120 kg
Weight: 14 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

**BARIATRIC. EASY WALKER. FOLDING WITH
FOREARM SUPPORTS. CABLE BRAKES. 4 WHEEL.
WM1018B**

Lightweight Aluminium
Adjustable height padded forearm supports
Adjustable length handbrakes
Total Width: 67 cm
Total Depth: 97 cm
Wheel Size: 20 cm (8")
Handle Height: 102 – 125 cm
Max. Weight Capacity: 180 kg
Weight: 10 kg



**TUTOR WALKER
CAN BE SUPPLIED WITH GUTTER ARMS OR
PADDED CHEST REST
WM1017TW**

Various Sizes
Two Models: 4 Wheel or
2 Wheel & Skis

More information upon request



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

2 WHEEL FRAMES / PICK-UP FRAMES

LIGHTWEIGHT FOLDING PICK-UP FRAME

WM1006

4 x Ferrules
Total Width: 59 cm
Handle Height: 80-91 cm
Width Between Handles: 43 cm
Max. Weight Capacity: 100 kg
Weight: 2.6 kg



BARIATRIC FOLDING FRAME. DUAL FRONT WHEELS (4) WITH SKIS / FERRULES

WM1006B

Dual front 5" wheels
Rear Skis / Ferrules
Total Width: 77.5 cm
Handle Height: 84 – 110 cm
Width Between Handles: 50 cm
Max. Weight Capacity: 250 kg
Weight: 4.5 kg



BARIATRIC FOLDING PICK-UP FRAME

WM1006BPU

4 x Ferrules
Total Width: 72 cm
Handle Height: 79 – 105 cm
Width Between Handles: 50 cm
Max. Weight Capacity: 250 kg
Weight: 3.7 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

RIGID NON-FOLDING ALUMINIUM PICK-UP FRAME **WM1007**

Total Width: 59 cm
Handle Height: 79 – 97 cm
Width Between Handles: 43 cm
Max. Weight Capacity: 113 kg
Weight: 2.2 kg



ULTRA-LIGHTWEIGHT, FOLDING, 2 WHEEL WALKING FRAME WITH SKIS / FERRULES **WM1009**

4 Sizes
The ultimate in simple folding without the use of buttons, levers, catches or other mechanisms
Light and easy to store and transport
Wheel Size: 12.5 cm (5")
Max. Weight Capacity: 125 kg
Total Weight: 2 kg



WM1009SUL
WM1009S
WM1009M
WM1009T

Handle Height
66 – 73 cm
74 – 81 cm
80 – 88 cm
85 – 93 cm

Total Width
65 cm
65 cm
66 cm
67 cm



Can easily be converted to a Pick-Up frame



SKIS - 2 STYLES

WM1904W1 - 22 mm (7/8") - Internal fit
WM1904W2 - 25 mm (1") - Internal fit
WM1904G1 - 25 mm (1") - External fit
WM1904G2 - 28 mm - External fit



FERRULES

Various Sizes: 13mm, 16mm, 19mm, 22mm, 25mm & 28mm



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

KNEE WALKER

WM1015

Fold-Down Steering Column
Adjustable Knee Height: 43 – 52 cm
Handle Height: 78 – 92 cm
Wheel Size: 20 cm
Max. Weight Capacity: 130 kg
Weight: 12 kg



TRAY WALKER – CABLE BRAKES

WM1014

Features 2 trays
Ideal for indoor use
Total Width: 55 cm
Total Depth: 55 cm
Wheel Size: 15 cm (6")
Handle Height: 76 – 92 cm
Max. Weight Capacity: 125 kg
Weight: 8.35 kg



LET'S GO INDOOR WALKER

WM1016

Folding, handle height adjustable, kitchen trolley for use indoors
Aluminium and wood frame with soft rubber handle
Removable plastic tray and vinyl bag
Can be operated one-handed
Swivel castors enable trolley to be turned in narrow spaces
Total Width: 56 cm
Total Depth: 70 cm
Wheel Size: 15 cm (6")
Tray Height: 70 cm
Handle Height: 85 – 95 cm
Weight: 7 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

ROOMBA INDOOR WALKER

WM1016D

Indoor Rollator/Walker

Walking frame is designed in contrasting white and black and has been developed for indoor use.

Features:

- Designed in contrasting black and white and developed for indoor use
- Full length handlebar is made from soft foamed rubber to provide a comfortable grip
- Full length brake handle suitable for one hand operation
- Very narrow with a width of only 57cm (22") ideally suited for use in tight spaces
- Foldable for easy storage and transportation and locks when folded up for easy handling
- Lightweight aluminium frame
- Includes a practical tray with a built-in cup holder to avoid spillages, and raised edges to prevent dropping items
- Includes basket to keep essentials handy

- Total Width: 57 cm
- Total Depth: 55 cm
- Height Min: 81cm
- Height Max: 94cm
- Weight: 6kg
- Weight Capacity: 100kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

UNDERARM ALUMINIUM CRUTCHES

		User Height
WM1060	TALL	178 - 198 cm
WM1061	MEDIUM	157 - 178 cm
WM1065	SMALL	137 - 157 cm
WM1072	CHILD	120 - 135 cm
WM1072S	PAEDIATRIC	106 - 131 cm

Superior quality lightweight anodised aluminium
Reinforcement under centre tube provides more load-bearing strength where it is most needed
Non-slip rubber tips provide good traction and built-in metal ring prevents excessive wear and ensures maximum durability
Max. Weight Capacity: 150 kg
Weight: 1.8 kg



ALSO AVAILABLE IN BARIATRIC MODELS

FOREARM ALUMINIUM CRUTCHES

		Sizing (Handle to Floor)
WM1062	TALL	84 - 107 cm
WM1063	MEDIUM	71 - 94 cm
WM1066	SMALL	56 - 79 cm
WM1073	PAEDIATRIC	47 - 70 cm
WM1073S	MINI PAEDIATRIC	40 - 55 cm

Pivoting cuff bracket is secured to a composite post for added durability
Standard handgrips, or ergonomic handgrips which help to distribute pressure on palms
Leg & forearm sections adjust independently for a custom fit
Height adjustable
Max. Weight Capacity: 110 - 150 kg
Weight: 1.9 kg



ALSO AVAILABLE IN HEIGHT ADJUSTABLE BARIATRIC MODEL STAINLESS STEEL - 200KG

Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

WALKING STICKS

CROOK HANDLE – ALUMINIUM. ADJUSTABLE **WM1048**

Handle Height: 76 – 99 cm
Max. Weight Capacity: 113 kg



OFFSET HANDLE – ALUMINIUM. ADJUSTABLE **WM1056**

Contoured handle with handgrips provide a comfortable and secure grip with convenient wrist strap
Attractive finish resists marring and scratching
Secured with a locking ring for added safety
Handle height adjustment in 2.5 cm increments: 76 – 99 cm
Max. Weight Capacity: 113 kg



BARIATRIC OFFSET HANDLE – STEEL. ADJUSTABLE **WM1056B**

Handle Height: 74 – 97 cm
Max. Weight Capacity: 225 kg
Colour: Black

* Same design as WM1056

T-SHAPE HANDLE - ALUMINIUM. ADJUSTABLE **WM1057BRZ – BRONZE**

Attractive finish resists marring and scratching
Push-button adjustment, circular ring marked on telescoping tube to indicate location of push-button
Secured with locking ring for added safety
Handle height adjustment in 2.5 cm increments: 75 – 98 cm
Non-slip rubber tips with metal inserts ensure a secure grip and a long life



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

WM1057BRZNS – BRONZE. NON-SLIP YELLOW TOP HANDLE

Handle Height Adjustment: 75 – 98 cm

ALSO IN ARTHRITIC GRIP STYLE



FOLDING – ALUMINIUM. ADJUSTABLE

WM1070A - SHORT

Handle Height Adjustment: 65 - 75cm (26"-30")



WM1070B - LONG

Handle Height Adjustment: 82 – 92cm (33"-37")

4 Folding sections fold out automatically
Adjusts in 2.5 cm increments



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

QUAD STICKS – ALUMINIUM. ADJUSTABLE

WM1054 - SMALL BASE (76 - 98 cm)

WM1055 - LARGE BASE (75 - 98 cm)

Sturdy aluminium lightweight construction

Push-button height adjustment with locking ring allow quick secure fitting

Non-slip rubber tips with metal insert ensure a secure grip and long life

Attractive finish resists marring and scratching

Handle height adjustment in 2.5 cm increments: 75 – 98 cm

Weight: 0.9 kg

WM1054



WM1055



PYRAMID STICK

WM1055P

Handle Height Adjustment: 75 – 95 cm

WM1055P



BARIATRIC QUAD – STEEL. SMALL BASE

WM1054B

BARIATRIC QUAD – STEEL. LARGE BASE

WM1055B

Handle Height Adjustment: 75 – 98 cm for both models

Max. Weight Capacity: 225 kg

Colour: Black

* Same designs as 1054 and 1055

Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

WHEELCHAIRS

STEEL - SELF-PROPELLING. PNEUMATIC TYRES

WM2020A – 46 CM

Developed for hospital and nursing home environment

Reversible, removable armrest

Allows for either a desk-arm or long armrest positioning serving the purpose of being able to be wheeled in to a dining table

(desk setting) or for easy access into or out of the wheelchair (long arm setting)

Removable, swing-away footplates

Spoke wheels

Heavy duty vinyl upholstery

Seat Width: 46 cm

Seat Depth: 40 cm

Seat Height: 50 cm

Back Height: 40 cm

Total Width: 68 cm

Wheel Size: 24" x 1 3/8"

Castor Size: 20 cm

Max. Weight Capacity: 130 kg

Weight: 19.6 kg



STEEL - SELF-PROPELLING. SOLID PU TYRES

SWINGAWAY, REMOVABLE, FOOTPLATES. REMOVABLE ARMRESTS

WM2020B – 46 CM

Seat Width: 46 cm

Seat Depth: 40 cm

Seat Height: 50 cm

Back Height: 40 cm

Total Width: 68 cm

Wheel Size: 24" x 1 3/8"

Castor Size: 20 cm

Max. Weight Capacity: 130 kg

Weight: 19.5 kg

* Same design as WM2020A except for tyre type

ALSO BARIATRIC MODELS

STEEL – SELF-PROPELLING. SOLID PU TYRES

SWINGAWAY, REMOVABLE, FOOTPLATES. REMOVABLE ARMRESTS

WM2040 - 50 CM

Seat Width: 50 cm

Seat Depth: 46 cm

Seat Height: 50 cm

Back Height: 40 cm

Total Width: 74 cm

Wheel Size: 24" x 1 3/8"

Castor Size: 20 cm

Max. Weight Capacity: 180 kg

Weight: 26.3 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

**STEEL – SELF-PROPELLING. SOLID PU TYRES
SWINGAWAY, REMOVABLE, FOOTPLATES. REMOVABLE
ARMRESTS**

WM2047 – 60 CM

Seat Width: 60 cm
Seat Depth: 46 cm
Seat Height: 50 cm
Back Height: 40 cm
Total Width: 82 cm
Wheel Size: 24" x 1 3/8"
Castor Size: 20 cm
Max. Weight Capacity: 200 kg
Weight: 24.3 kg



STEEL TRANSIT WHEELCHAIR. SOLID PU TYRES

WM2030 – 46 CM

Removable / Swing away footplates
Removable padded armrests
Heavy duty upholstery
Seat Width: 46 cm
Seat Depth: 40 cm
Seat Height: 50 cm
Back Height: 40 cm
Total Width: 60 cm
Wheel Size: 12 1/2" x 2 1/4"
Castor Size: 20 cm
Max. Weight Capacity: 130 kg
Weight: 17.5 kg



**STEEL TRANSIT WHEELCHAIR. BARIATRIC. SOLID PU
TYRES. SWINGAWAY, REMOVABLE, FOOTPLATES.
REMOVABLE ARMRESTS**

WM2040T – 50 CM

Seat Width: 50 cm
Seat Depth: 40 cm
Seat Height: 50 cm
Back Height: 40 cm
Total Width: 64 cm
Wheel Size: 12 1/2" x 2 1/4"
Castor Size: 20 cm
Max. Weight Capacity: 180 kg
Weight: 18.5 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

ACCESSORIES FOR STEEL WHEELCHAIRS



IV POLE
WMN2058



ELEVATING LEG RESTS

WM1906EL – SELF PROPELLING
WM1909EL – TRANSIT MODEL

WHEELCHAIR OXYGEN CYLINDER HOLDER (FOR STEEL WHEELCHAIRS)

WM2056W

Attaches easily to wheelchair
Also model to suit Karma wheelchairs



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

KARMA EAGLE LIGHTWEIGHT WHEELCHAIR

EAGLE 14" x 16"

WM2006-36

Seat Width: 36 cm
Seat Depth: 41 cm
Seat Height: 50 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.2 kg

EAGLE 16" x 16"

WM2006-40

Seat Width: 41 cm
Seat Depth: 41 cm
Seat Height: 50 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.4 kg

EAGLE 17" x 16"

WM2006-43

Seat Width: 43 cm
Seat Depth: 41 cm
Seat Height: 50 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.5 kg

EAGLE 18" x 16"

WM2006-46 (x40)

Seat Width: 46 cm
Seat Depth: 41 cm
Seat Height: 50 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.6 kg

EAGLE 18" x 18"

WM2006-46 (x46)

Seat Width: 46 cm
Seat Depth: 46 cm
Seat Height: 50 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.6 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

EAGLE 19" x 18"**WM2006-49**

Seat Width: 49 cm
Seat Depth: 46 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 160 kg
Weight: 15.5 kg

EAGLE 20" x 18"**WM2006-50**

Seat Width: 52 cm
Seat Depth: 46 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 160 kg
Weight: 15.7 kg

EAGLE 22" x 18"**WM2006-55**

Seat Width: 55 cm
Seat Depth: 46 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 160 kg
Weight: 15.9 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

KARMA FLEXX SELF-PROPELLING WHEELCHAIR
WM2001F-46

Seat Width: 46 cm
Seat Depth: 41 cm
Seat Height: 50 cm
Back Height: 43 cm
Max. Weight Capacity: 130 kg
Weight: 14.6 kg

SIZES AVAILABLE: (16" to 24")



KARMA TRANSIT II LIGHTWEIGHT WHEELCHAIR

TRANSIT II 16" x 16"
WM2008-40

Seat Width: 41 cm
Seat Depth: 41 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 13.5 kg

TRANSIT II 17" x 16"
WM2008-43

Seat Width: 43 cm
Seat Depth: 41 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 13.9 kg

TRANSIT II 18" x 16"
WM2008-46 (x40)

Seat Width: 46 cm
Seat Depth: 41 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.3 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

TRANSIT II 18" x 18"**WM2008-46 (x46)**

Seat Width: 46 cm
Seat Depth: 46 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 130 kg
Weight: 14.3 kg

TRANSIT II 20" x 18"**WM2008-50**

Seat Width: 51 cm
Seat Depth: 46 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 160 kg
Weight: 15.2 kg

**TRANSIT II 22" x 18"****WM2008-55**

Seat Width: 55 cm
Seat Depth: 46 cm
Seat Height: 52 cm
Back Height: 40 cm
Max. Weight Capacity: 160 kg
Weight: 15.7 kg

KARMA ERGO LITE DELUXE SELF – PROPELLING WHEELCHAIR**ERGO LITE DELUXE 16" x 16"****WM2011-40**

Seat Width: 41 cm
Seat Depth: 42 cm
Seat Height: 48 cm
Back Height: 42 cm
Max. Weight Capacity: 100 kg
Weight: 10.2 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

Page 24

ERGO LITE DELUXE 18" x 16"
WM2011-46

Seat Width: 45 cm
Seat Depth: 42 cm
Seat Height: 48 cm
Back Height: 42 cm
Max. Weight Capacity: 100 kg
Weight: 10.3 kg



KARMA S-ERGO 125 SELF-PROPELLING WHEELCHAIR

S-ERGO 125 16" x 17"
WM2000-40

Seat Width: 41 cm
Seat Depth: 43 cm
Seat Height: 51 cm
Back Height: 40 cm
Max. Weight Capacity: 115 kg
Weight: 14.2 kg

S-ERGO 125 18" x 17"
WM2000-46

Seat Width: 45 cm
Seat Depth: 43 cm
Seat Height: 51 cm
Back Height: 40 cm
Max. Weight Capacity: 115 kg
Weight: 14.5 kg

S-ERGO 125 20" x 17"
WM2000-50

Seat Width: 51 cm
Seat Depth: 43 cm
Seat Height: 51 cm
Back Height: 40 cm
Max. Weight Capacity: 115 kg
Weight: 14.8 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

KARMA ERGO LITE TRANSIT WHEELCHAIR

ERGO LITE 16" x 16"

WM2007-40

Seat Width: 41 cm
Seat Depth: 42 cm
Seat Height: 48 cm
Back Height: 42 cm
Max. Weight Capacity: 100 kg
Weight: 8.1 kg

ERGO LITE 18" x 16"

WM2007-46

Seat Width: 45 cm
Seat Depth: 42 cm
Seat Height: 48 cm
Back Height: 42 cm
Max. Weight Capacity: 100 kg
Weight: 8.3 kg



KARMA ERGO LITE DELUXE TRANSIT WHEELCHAIR

ERGO LITE DELUXE 16" x 16"

WM2007D-40

Seat Width: 41 cm
Seat Depth: 42 cm
Seat Height: 48 cm
Back Height: 42 cm
Max. Weight Capacity: 100 kg
Weight: 8.7 kg

ERGO LITE DELUXE 18" x 16"

WM2007D-46

Seat Width: 45 cm
Seat Depth: 42 cm
Seat Height: 48 cm
Back Height: 42 cm
Max. Weight Capacity: 100 kg
Weight: 8.8 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

KARMA S-ERGO 125 TRANSIT WHEELCHAIR

S-ERGO TRANSIT 16" x 17"

WM2000T-40

Seat Width: 41 cm
Seat Depth: 43 cm
Seat Height: 51 cm
Back Height: 40 cm
Max. Weight Capacity: 115 kg
Weight: 13.1 kg



S-ERGO TRANSIT 18" x 17"

WM2000T-46

Seat Width: 45 cm
Seat Depth: 43 cm
Seat Height: 51 cm
Back Height: 40 cm
Max. Weight Capacity: 115 kg
Weight: 13.1 kg

S-ERGO TRANSIT 20" x 17"

WM2000T-50

Seat Width: 51 cm
Seat Depth: 43 cm
Seat Height: 51 cm
Back Height: 40 cm
Max. Weight Capacity: 115 kg
Weight: 13.3 kg

KARMA TILT WHEELCHAIR

KM-VIP-515 16" x 17"

WM2012-40

Seat Width: 41 cm
Seat Depth: 43 cm
Seat Height: 50 cm
Back Height: 41 cm
Max. Weight Capacity: 115 kg
Weight: 17.5 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

KM-VIP-515 18" x 17"**WM2012-46**

Seat Width: 45 cm
Seat Depth: 43 cm
Seat Height: 50 cm
Back Height: 41 cm
Max. Weight Capacity: 115 kg
Weight: 17.5 kg

KARMA RECLINE**WM2013**

Seat Width: 41 cm
Seat Depth: 41 cm
Seat Height: 51 cm
Back Height: 48 cm
Max. Weight Capacity: 115 kg
Weight: 16 x 16" 18.3 kg
18 x 18" 19 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

KARMA ACCESSORIES



IV POLE



OXYGEN HOLDER



BACKREST BAG



ANTI-TIPPER



DINNER TRAY



ELEVATING FOOTREST



STUMP SUPPORT



FOLDABLE HEADREST

SOMA ECON LIGHTWEIGHT SELF-PROPELLING WHEELCHAIR. SOLID PU TYRES

Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

WM2009-40

Seat Width: 41 cm
Seat Depth: 43 cm
Seat Height: 48 cm
Back Height: 43 cm
Total Width: 59 cm
Wheel Size: 24" x 1 3/8"
Castor Size: 20 cm
Max. Weight Capacity: 115 kg
Weight: 14.6 kg

**SOMA ECON LIGHTWEIGHT SELF-PROPELLING WHEELCHAIR.
SOLID PU TYRES****WM2010-46**

Seat Width: 45 cm
Seat Depth: 43 cm
Seat Height: 48 cm
Back Height: 43 cm
Total Width: 64 cm
Wheel Size: 24" x 1 3/8"
Castor Size: 20 cm
Max. Weight Capacity: 115kg
Weight: 14.8 kg

* Same design as 2009

**SOMA ECON LIGHTWEIGHT TRANSIT WHEELCHAIR. SOLID PU TYRES
WM2031-40**

Seat Width: 41 cm
Seat Depth: 43 cm
Seat Height: 48 cm
Back Height: 43 cm
Total Width: 55.5 cm
Wheel Size: 16" x 1.75"
Castor Size: 20 cm
Max. Weight Capacity: 115 kg
Weight: 14.1 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

SOMA ECON LIGHTWEIGHT TRANSIT WHEELCHAIR. SOLID PU TYRES
WM2032-46

Seat Width: 45 cm
Seat Depth: 43 cm
Seat Height: 48 cm
Back Height: 43 cm
Wheel Size: 16" x 1.75"
Castor Size: 20 cm
Max. Weight Capacity: 115 kg
Weight: 14.3 kg

* Same design as 2031

SPARE PARTS AVAILABLE: TYRES, TUBES ETC.

TRANSFER BOARDS

Assorted Range



RAMPS

PORTABLE ROLL-UP WHEELCHAIR RAMP

WM2050

Can be easily dismantled and rolled up
Length: 90 cm
Width: 75 cm
Max. Weight Capacity: 250 kg
Weight: 7 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

PRIDE MODEL – SILVER STAR. ALUMINIUM
WM2050-1220 1.2 Metre

Width: 75 cm
Length: 122 cm
Max. Height: 30 cm
Max. Weight Capacity: 272 kg
Weight: 11 kg



PRIDE MODEL – SILVER STAR. ALUMINIUM
WM2050-1830 1.8 Metre

Width: 75 cm
Length: 183 cm
Max. Height: 51 cm
Max. Weight Capacity: 272 kg
Weight: 16 kg

DECPAC - FIBREGLASS WITH EDGE BARRIER LIMITER – 1.2 METRE
WM2051

Width: 82 cm
Length: 120 cm
Max. Height: 30 cm
Max. Weight Capacity: 300 kg
Weight: 9.7 kg

Full range including EBL 900–3500 MM



RAMP / WEDGE (RUBBER)
WM2050RW

100% recycled rubber, will not slip, reducing the risk of trips and falls
The tough, long lasting rubber creates a smooth transition for wheelchairs and for those using walking frames
Made to order Price dependent on size



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au


WHEELCHAIR CUSHIONS

Various sizes and weight capacities including Equagel, Roho, Jay, Diffuser, Multi-Purpose

EQUAGEL GENERAL – 3.2 CM THICKNESS

WM5008

40 x 40 cm (16 x 16")
46 x 40 cm (18 x 16")
46 x 46 cm (18 x 18")
50 x 46 cm (20 x 18")
Also wider sizes available up to 76 x 50 cm (30 x 20")



EQUAGEL PROTECTOR – 5 CM THICKNESS

WM5009

40 x 40 cm (16 x 16")
46 x 40 cm (18 x 16")
46 x 46 cm (18 x 18")
50 x 40 cm (20 x 16")
50 x 46 cm (20 x 18")
Also wider sizes available up to 76 x 50 cm (30 x 20")



EQUAGEL ADJUSTABLE PROTECTOR – 6.4 CM

WM5010

Three different firmnesses of 3.2 cm EquaGel are provided for the purpose of selecting the perfect combination

40 x 40 cm (16 x 16")
46 x 40 cm (18 x 16")
46 x 46 cm (18 x 18")
50 x 40 cm (20 x 16")
50 x 46 cm (20 x 18")
Also wider sizes available up to 76 x 50 cm (30 x 20")

ROHO CUSHIONS



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

JAY CUSHIONS



DIFFUSER CUSHION. MEMORY FOAM **WM5000**

Medium: 44 x 42 cm
Large: 49 x 46 cm



MULTI-PURPOSE CUSHION

Unique egg foam surface diffuses shearing forces.
Also fits wheelchairs. Premium-grade foam ensures product longevity.
Ideal for sitting for long periods.



KEY HOLE CUSHION

Encourages correct pelvic tilt to reduce pressure on the lumbar discs and thus minimise potential back problems.
Symmetrical aperture means tender parts are left 'free' while the 9-degree angle helps maintain your natural lumbar curve.
Designed to take the weight on the user's body contact points (ischium) and along thighs.



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

RING RELIEF CUSHION

For post-natal and post-operative discomfort, haemorrhoids, bruised coccyx and other painful afflictions.
Ideal compound to diffuse and spread the pressure.
Will not sag or lose its shape.
Covered in moisture-resistant bacteria inhibiting (hospital grade) Steri-Plus.



POSTURE WEDGE

Promotes better posture – helps maintain your normal lumbar curve while sitting.
Relieves stress on back, increases user comfort for longer periods.



LEG RELAXER

Recommended for those with varicose veins, swollen feet and ankles, 'tired' legs, during pregnancy and convalescing.
Even weight distribution alleviates lower back strain.
Minimises tension in the sciatic nerve.
Adjustable for two levels of comfort.



BED WEDGE

Mildly angled to minimise pressure to the weight bearing lower vertebrae.
Encourages correct spinal alignment spreads the load, provides body-hugging support.
Ideal for persons suffering medical conditions requiring an angled sleeping position.



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

BACK SUPPORTS

LUMBAR ROLLS

Available in various shapes.
Helps prevent the onset of lower back pain.

D-Shape: Regular size, form-fits car seats, office chairs and couches.

14 x 9 x 31 cm long.

Round: Circular. Compresses to form ellipse.

11 cm diameter x 30 cm long.

Oval: Unique rippled surface diffuses pressure points.

17 x 7 x 35 cm long.

Egg Shape: Accommodates medium to extreme lower back curvature.

13 x 15 x 34 cm long.



CONTOURED FULL BACK SUPPORT

Ideal for persons with middle and lower-back conditions.
Therapeutically designed to comfort and support the spine.
Helps maintain the spine's natural curvature.



FLO BAC ADJUSTABLE CAR – BACK / SEAT SUPPORT WM5003

Compensates for poor seating – less stress on your back,
less fatigue.

Lightweight and portable.

Durable open-weave soft mesh cover allows air to flow
– 'no sweat'.

Adjustable mesh with 'massaging' beaded wood seat attachment.



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

KOALA COMFORT BACK SUPPORT

Adjustable internal bladder adapts to the individual & allows stressless alignment of your back.

Supportive flanks for lateral stability.

Absorbs body pressures – helps alleviate muscle fatigues.



TOTAL SPINAL SUPPORT

Ensures your spine is fully supported at all times.

Ideal for persons with middle and lower-back conditions.

Contoured to fit snugly into your back.

Also features flanks which inhibit lateral movement.

Stable on any seat or chair.

Model 1: 40 x 55 cm high – For use on firm seating (car, office chair)

Model 2: 40 x 50 cm high – For use on soft seating (couch / home)



SHEEPSKINS / LAMBSWOOL

Heel / Elbow Protectors

Overlay

Wheelchair Cover

Etc



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

SCOOTERS

CELEBRITY DX **WM8000**

High performance scooter – 4 Wheel
Ideal for outdoor use
Front-mounted high back swivel seat with headrest and armrests
2 x 12V Rechargeable batteries
Pneumatic tyres
Easily adjusted steering mount with LED headlight and indicator lights
Front basket
Length: 119 cm
Width: 61 cm
Max. Speed: Up to 9.5 km/h
Range: 40 km
Ground Clearance: 8.2 cm
Turning Radius: 150 cm
Max. Weight Capacity: 159 kg
Total Weight: 97 kg (includes batteries – 22 kg)



AND MORE UPON REQUEST

GOGO ULTRA X **WM8001**

Available in 3 or 4 wheel configuration
Disassembles into 5 pieces for easy transportation and storage
Front-mounted post seat with armrests
Easy to remove rechargeable battery
Adjustable tiller
Flat free tyres
Front basket
Length: 101 cm (4 Wheel) 94 cm (3 Wheel)
Width: 49.5 cm (both models)
Max. Speed: Up to 6.4 km/h
Range: 11 km
Ground Clearance: 3.2 cm (both models)
Turning Radius: 112 cm (4 Wheel) 84cm (3 Wheel)
Max. Weight Capacity: 118 kg
Weight (without batteries): 36.3 kg (4 Wheel) 34.5 kg (3 Wheel)



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

BEDROOM

SELF-HELP POLE – FREESTANDING

WM3012

The free-standing gooseneck is fabricated from heavy tubular steel section and is epoxy powder coated

Chain link adjustment gives a handle height variation of approximately 500 mm

Max. Weight Capacity: 120 kg



BED POLE – CROOK HANDLE

WM4000CR

Fits under mattress & provides firm handle to assist with moving & positioning in bed

Suitable for use on single, double, queen & king size beds

Pole Height: 63 cm



DELUXE ALUMINIUM I.V. STAND

WM4001

Circular "Ram Horn" style hooks

Corrosion-resistant aluminium stem with heavyweight steel base

Friction-grip lock permits height adjustments without danger of slipping

Excellent stability with five-legged enamelled steel base

Height adjusts from 144 – 243 cm



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

Page 39

BED CRADLE

WM4002

Sturdy steel construction
Supports bed linen to keep pressure off patient's legs and feet
Easily fits under mattress on either end or side of bed



SINGLE BED STICK – DOUBLE PRONG

WM4003

Fits under mattress and provides firm handle to assist with moving and positioning in bed
Can be used as a safety rail



CHAIRS

SMALL BACK SUPPORT CHAIR

WM6000

Fully upholstered support chair
Sculptured backrest for maximum support and comfort
Powder coated steel frame
Easy to clean vinyl
Seat Width: 46 cm
Seat Depth: 51 cm
Seat Height: 48 -60 cm
Back Height: 35 cm
Total Width: 63 cm
Max. Weight Capacity: 120 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

Page 40

BARIATRIC SMALL BACK SUPPORT CHAIR

WM6000B

Seat Width: 56 cm
Seat Depth: 50cm
Seat Height: 50 – 61 cm
Back Height: 34 cm
Total Width: 73 cm
Max. Weight Capacity: 210 kg

* Same design as WM6000

MID-BACK SUPPORT CHAIR

WM6000MB

Mid-back chair offering greater lumbar support due to increased backrest height
Seat Width: 45 cm
Seat Depth: 50 cm
Seat Height: 48 - 60 cm
Back Height: 55 cm
Total Width: 63 cm
Max. Weight Capacity: 120 kg



COMFORT CHAIR

WM6001HB

Multi-purpose chair featuring a large contoured backrest for both lumbar and neck support
Seat Width: 57 cm
Seat Depth: 47 cm
Seat Height: 48 - 60 cm
Back Height: 67 cm
Total Width: 73 cm
Max. Weight Capacity: 120 kg

Also available in bariatric size: WM6001B (210 kg)



FOOTSTOOL

WM6002

Angled, height adjustable footstool (Optional extra for WM6001HB)
Height Adjustment: 43 - 55 cm

Other models available



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

TODD CHAIR
WM6004

Seat height, seat length and neck rest are easily adjusted
The side build up, inbuilt lumbar support and seat bolster ensure correct postural support
The wide padded armrests give comfortable shoulder support and ease of rising
Seat Width: 52 cm
Seat Depth: 41 - 51 cm
Seat Height: 51 – 58 cm
Back Height: 75 cm
Total Width: 67 cm
Max. Weight Capacity: 135 kg
Mobility castors fitted



OVERBED TABLE
WM6052

"C" base positions around chair or under bed enabling close positioning to the user
Spring loaded height adjustment for ease of use for carers and users
Various sizes and designs available



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

CONTINENCE PRODUCTS

Bed Pads
Linen Savers
Chair Pads
Wheelchair Pads



LIFTING HOISTS / SLINGS

All popular models, including standing lifters, and full range of slings can be supplied upon request



BEDSIDE COMMODE - ECONOMY WM3000

Made from a mild steel, epoxy powder coat
Phosphate pre-treated frame
Lift up toilet seat and fixed bowl carrier
Padded vinyl back and removable padded seat
Seat Width: 44 cm
Seat Depth: 44 cm
Seat Height: 49 – 61cm
Back Height: 44 cm
Total Width: 66 cm
Max. Weight Capacity: 110 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

BEDSIDE COMMUNE – DELUXE

WM3007

Fully upholstered adjustable height commune
Designed to be discreet and stylish
Complete with hinged seat, back, arms and skirt to conceal the commune pan
Seat Width: 45 cm
Seat Depth: 50 cm
Seat Height: 39 – 52 cm
Total Width: 63 cm
Max. Weight Capacity: 110 kg
Weight: 12.95 kg



BARIATRIC MODELS AVAILABLE 175 - 210 KG

COMMUNE / OVERTOILET FRAME

WM3001

Aluminium
Full length plastic armrests
Seat Width: 46 cm
Seat Depth: 38 cm
Seat Height: 40 – 50 cm
Total Width: 66 cm
Max. Weight Capacity: 136 kg
Weight: 3.9 kg



WEIGH SCALE CHAIR – DIGITAL

WM3998

Dry battery or AC adaptor power
Comfortable, white plastic seat
Rugged and reliable operation
Full digital calibration ensuring accurate and simple calibration
Auto power off to preserve battery life
Safety wheels with brakes
Max. Weight Capacity: 150 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

PRINCESS AIR CHAIR / DELUXE AIR BED

WM6003

Utilising the unique air comfort seating system
Complete workstation for the nursing of dependant patients offers proven pressure relief
Total surface air pressure system
No gaps in seating surface
Adjustable seat (infinite locking)
Adjustable back (infinite locking)
Adjustable wings and drop-down arms
Adjustable leg rest (infinite locking)
Detachable / adjustable head rest
Locking castors (standard and directional)
Full metal frame
Replaceable cover

Various models



HOSPITAL BED - VARIOUS MODELS

MATTRESSES

COMMUNITY CARE BED

WM3050

Community electric care bed. 4 Section - pull apart
Electronically controlled bed with self-help pole for aged care or home environments
Enables both care giver and patient to reposition the bed for most aspects of care
Four function – electric high/low, backrest, knee bend, Trendelenburg & reverse Trendelenburg
External Length: 2200 mm
External Width: 930 mm
Bed Height: 340 – 750 mm
Mattress dimensions: 2000 x 850 mm
Max. Weight Capacity: 120 kg




LIFE COMFORT MATTRESS

WM3051

Dimensions: 2000 x 850 x 130 mm

Many other sizes



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

BARIATRIC KING SINGLE BED
WM3050KS

Bariatric 2 or 3 function electric lifting and reclining bed
External Length: 2120 mm
External Width: 1070 mm
Bed Height Range: 200 – 750 mm
Mattress Dimensions: 2030 x 1050 mm
Max. Weight Capacity: 350 kg



POWERLIFT RECLINER CHAIRS

PRIDE C5
WM6010

3-Position: Recline, Lift & Leg Raise
Seat to Floor: 50 cm
Depth: 51 cm
Width: 56 cm
Back Width: 71 cm
Backrest Height: 63 cm
Max. Weight Capacity: 170 kg
Integrated Battery Backup: Yes
Colour: Arctic Blue



PRIDE C101
WM6011

3-Position: Recline, Lift & Leg Raise
Seat to Floor: 47 cm
Depth: 51 cm
Width: 51 cm
Back Width: 72 cm
Backrest Height: 70 cm
Max. Weight Capacity: 147 kg
Integrated Battery Backup: Yes
Colour: Sandal

Other models available



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

BATHROOM

OVER TOILET FRAME

WM3008

Powder coated corrosion resistant galvanised tubing frame
Deluxe commode pan enables the product to also function as a commode

Seat Width: 45 cm

Seat Depth: 43 cm

Seat Height: 48 – 60 cm

Max. Weight Capacity: 110 kg

Weight: 3.2 kg

Various Sizes: (200 kg)

WM3011 FOLDING MODEL ALSO AVAILABLE



BARIATRIC OVER TOILET FRAME

WM3008B

Seat Width: 53 – 57 cm tapered

Seat Height: 42.5 – 55.5 cm

Max. Weight Capacity: 200kg

Weight: 4.25 kg



BARIATRIC COMMODE / OVER TOILET FRAME

WM3003B

All steel frame with Versa Guard coating is durable and warm to touch

Comes with bucket and lid

Easy to operate drop arm mechanism facilitates lateral transfer

Seat Width: 63 cm

Seat Depth: 46 cm

Seat Height: 50 – 60 cm

Max. Weight Capacity: 225 kg

Weight: 9.85 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

SHOWER CHAIR**WM3004**

Powder coated corrosion resistant galvanised tubing frame
Plastic moulded seat allowing for both cleanliness and comfort
Seat Width: 45 cm
Seat Height: 42 – 54 cm
Total Width: 51 cm
Max. Weight Capacity: 110 kg
Weight: 3.2 kg
Other sizes available

**BARIATRIC ADJUSTABLE SHOWER CHAIR****WM3004B**

Seat Width: 45 cm
Seat Height: 44 – 57 cm
Total Width: 68 cm
Max. Weight Capacity: 200 kg
Weight: 4.55 kg

**SHOWER STOOL****WM3005**

Powder coated corrosion resistant galvanised tubing frame
Adjustable legs
Armrests assist with patient transfer on and off the chair ensuring safety and confidence
Other styles available including "corner" model (triangular shape)
See WM3005T
Seat Width: 47 cm
Seat Depth: 29 cm
Seat Height: 51 – 63 cm
Max. Weight Capacity: 110 kg
Weight: 2.9 kg



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

BARIATRIC SHOWER STOOL
WM3005B

Seat Width: 52 cm
Seat Depth: 30 cm
Seat Height: 47 – 59 cm
Max. Weight Capacity: 200 kg
Weight: 4.45 kg



TRIANGULAR SHOWER STOOL
WM3005T

Seat Width: 52 cm
Seat Depth: 37 cm
Seat Height: 44 -54 cm
Max. Weight Capacity: 160 kg
Weight: 2.2 kg



SWIVEL BATH SEAT
WM3006

Powder coated stainless steel
Plastic moulded seat
Width between arms 48 cm
Max. Weight Capacity: 100 kg
Weight: 6 kg

Height adjustment spacers available: 3" & 6"



FOLDING SHOWER STOOL
WM3100

Curved plastic seat with handles
Drain holes help to keep water off sitting area
Aluminium frame designed to provide both comfort and safety
Folding for easy transport
Seat Width: 47 cm
Seat Depth: 29 cm
Seat Height: 36.8 – 39.4 cm
Total Width: 47 cm
Total Depth: 54 cm
Max. Weight Capacity: 110 kg
Weight: 2.2 kg



BATH MAT

Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

WM3500

SHOWER MAT

WM3501

URINALS

WM3503

Spillproof, Male & Female



BED PANS

WM3504

Slipper / Standard



SHOWER COMMUNE TRANSPORTER

WM3105

Stainless steel construction
Drop-down arms
Locking castors
Seamless padded seat and backrest
Bed pan holder
Seat Width: 44 cm
Seat Depth: 43 cm
Seat Height: 60 cm
Width Between Arms: 45 cm
Max. Weight Capacity: 125 kg
Weight: 15.5 kg
Optional extras: Oxygen cylinder holder & I.V. Pole

Various models including self-propelling



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

SHOWER COMMODE TRANSPORTER - ATLANTIC WAVE
WM3107

High grade silver painted aluminium frame
Swing-up armrests
Polyurethane backrest and seat
Removable over seat
Fold-up footrest
4 Wheel brake castors
Seat Width: 43 cm
Seat Depth: 46 cm
Seat Height: 56 cm
Width Between Armrests: 45 cm
Max. Weight Capacity: 160 kg
Weight: 11.5 kg



BATH BENCH
WM3102

Non-sliding
4 Adjustable Legs
Backrest
Seat Width: 67 cm
Seat Depth: 41 cm
Total Width: 77 cm
Total Depth: 53 cm



BATH BOARD
WM3101

Provides seating while showering in the bath
Adjustable shower board with handgrip
Comfortable contoured seat for bathing
Anti-slip rubber pad improves safety and protects bath surfaces from scratching
Drain holes help to keep water off sitting area
Portable & easy to assemble
Width: 69 cm
Depth: 26 cm
Max. Weight Capacity: 150 kg
Weight: 3 kg
Other styles available



Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999
Email: enquiries@wheelchairman.com.au

RAISED TOILET SEAT

WM3002

Raised toilet seat with arms

Arm Height: 15 cm

Arm Width: 50 cm

Seat Raised: 12.5 cm



WM3002A

Seat Raised: 5 cm

WM3002B

Seat Raised: 10 cm

WM3002C

Seat Raised: 15 cm

WM3002A, WM3002B & WM3002C also available with lid



BATHTUB BAR

WM3200

White enamel-coated 2.5 cm steel bar provides reliable patient assistance into and out of bathtub

Crank style locking device attaches to bath walls up to 17 cm wide

Unit is 15 cm wide & extends to 35 cm above bath wall

Equipped with rubber pads to protect bath walls

Weight: 3 kg



REACHERS

Homecraft models

Various sizes



AIDS TO DAILY LIVING

Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

2 Handled Mug & Lids

Caring Mug

Handy Bars



Jar Opener



Key Turner



Rota Swivel Seat



PEDAL EXERCISER

WM7000

Smooth pedal movement
Great for exercising the knees
Good for circulation
Max. Weight Capacity: 113 kg
Weight: 2.5 kg



Long Handled Sponge

Plate Surround

Sock Aid

Stocking Aid

Toe Washer

This is just a small sample of the Daily Living products we supply.

If the product does not feature in this catalogue, please contact us at the Showroom (03 9312 2888) and we will source it for you.

Wheelchair Man (Australasia) Pty. Ltd.

Phone: (03) 9312 2888 Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

**473 Ballarat Road
Sunshine Victoria 3020**

Melways Reference. Map 26, J10



WHEELCHAIR MAN (AUSTRALASIA) PTY. LTD.

Jim Nolan - Director

ABN: 78 440 131 015

Phone: (03) 9312 2888

Fax: (03) 9312 2999

Email: enquiries@wheelchairman.com.au

Web: www.wheelchairman.com.au