



CATALOGUE

Walkers
Walking Frames
Walking Aids
Wheelchairs
Scoters
Bedroom
Powerlift Chairs
Bathroom
Aids to Daily Living



* Also Hire, Service & Repair of our Products

ADDRESS: 473 Ballarat Road

Sunshine Vic 3020

PHONE: (03) 9312 2888

FAX: (03) 9312 2999

EMAIL: enquiries@wheelchairman.com.au

WEB: www.wheelchairman.com.au

TABLE OF CONTENTS

Rollator 2

CONFIDENTIAL

ROLLATOR

AIRGO® FUSION™ SIDE-FOLDING ROLLATOR AND TRANSPORT CHAIR

SKU: WM1999

THE TWO IN ONE SIDE-FOLD ROLLING WALKER AND TRANSPORT CHAIR!

The Airgo Fusion was designed to allow you to remain independent and live your life to the fullest. Its lightweight design and unparalleled stability provide the confidence to walk safely again. And if you ever find that you've ventured out too far, the Airgo Fusion converts easily to a transport chair so that a loved one can help bring you back home in comfort and style. Find out more about why so many people love their Airgo Fusion!

Features:

- Lightweight side-folding rollator and transport chair in one design
- Folds sideways easily by just a pull of the strap!
- Patent pending Swing n' Lock footrests clip to the frame when folded
- The frame locks together when folded for safer lifting and storing
- Padded seat and a full backrest for superior comfort
- Dual armrests provide a place to rest yours arms in either mode
- Smooth bearings and an improved centre of gravity for easy pushing
- Unique ergonomic hand-grips fit naturally for easier brake access
- Features our new senior friendly tool-free brake adjustment knob
- Lightweight frame (17.1 lb net) designed for superior stability
- Reflective soft basket can be removed and is machine washable

Specifications:

- Seat width: 46 cm
- Seat Depth: 43 cm
- Total Height: 76 cm – 93 cm
- Seat to Floor: 51 cm
- Length: 72 cm, Width: 69 cm
- Suitable for a person: 153 - 188 cm in height
- Max. Weight Capacity: 136 kg
- Total Weight: 10.1 kg

Categories: [Specialist Walkers](#), [Standard Size](#), [Transit](#), [Walkers](#), [Walking Aids](#)





ASPIRE VOGUE, CARBON FIBRE, SEAT WALKER / ROLLATOR

SKU: WM1003VCA

ASPIRE VOGUE, CARBON FIBRE, SEAT WALKER / ROLLATOR

Features:

- The stylish European designed Aspire Vogue Walkers are a lightweight, versatile range with a sturdy frame and easy to fold locking system, suitable for compact transportation when folded
- Luxury Look and Feel - European design, profiled tubing, moulded fittings and chic colour palette
- Lockable Handbrakes - For increased user safety when walking and when sitting
- Lightweight and Ergonomic - Easy to push with a small turning circle and low rolling resistance
- Unbreakable TPE Super Soft Wheels - High durability with a smooth ride for increased comfort
- Locks when Folded - Makes lifting the frame easier and minimises storage space required
- Removable Utility Bag - Stylish and simple to remove for shopping

Specifications:

- Seat width: 46 cm
- Seat Depth: 34cm
- Total Height: 69.5cm - 96 cm
- Seat to Floor: 55 cm
- Length: 680 mm, Width: 600 mm
- Weight capacity: 150kg
- Total Weight: 5.7 kg

Categories: [Standard Size](#), [Walkers](#), [Walking Aids](#)

CONFIDENTIAL



WALKER. A-FRAME. STANDARD

SKU: WM1003MW

SPECIFICATIONS

- Seat Height: 55cm
- Seat width: 34.5cm
- Sitting space: 44cm
- Handle height: 80cm-91cm
- User height: 151-170cm
- Safe working load: 130kg

Categories: [Standard Size](#), [Walkers](#)

CONFIDENTIAL