



CATALOGUE

Walkers
Walking Frames
Walking Aids
Wheelchairs
Scoters
Bedroom
Powerlift Chairs
Bathroom
Aids to Daily Living



* Also Hire, Service & Repair of our Products

ADDRESS: 473 Ballarat Road

Sunshine Vic 3020

PHONE: (03) 9312 2888

FAX: (03) 9312 2999

EMAIL: enquiries@wheelchairman.com.au

WEB: www.wheelchairman.com.au

TABLE OF CONTENTS

tri-walker 2

CONFIDENTIAL

TRI-WALKER



WM1010D - TRI-WALKER INCLUDES BAG

SKU: WM1010D

Product Features

- Rustproof and robust aluminium frame ensures longevity and is ultra-light (only 4.7kg (10lb)) to manoeuvre / transport
- Handle height easily adjusts to suit various users
- 18cm (7") diameter solid, lightweight wheels move effortlessly indoor or outdoor and are low maintenance
- Weight Capacity: 115kg
- Front wheel swivels for increased manoeuvrability
- Arthritic friendly, looped cable brakes are adjustable, easily activated and lockable for safety
- Easy one hand folding for storage / transportation

Categories: [Specialist Walkers](#), [Walkers](#)

CONFIDENTIAL